

Product Name :
Say Bolt Viscometer Apparatus

Product Code :
LIM-CAT-L0001-00001



Description :

Say Bolt Viscometer Apparatus

Technical Specification :

The Saybolt Two-Tube Digital Viscometer is supplied complete with:

Universal Orifice, 2 pcs

Furol Orific, , 2 pcs

Heat Transfer Oil, 5 lt

Key

Saybolt viscosity flask, glass, 60 ml, 2 pcs.

Thermometer Support, 2 pcs

Stopper, 2 pcs

Technical Specifications :

Dimensions: 450x300x550 mm

Weight (approx.): 10 kg

Power: 750 W

Say Bolt Viscometer Apparatus, Say Bolt Viscometer Apparatus Bulk Suppliers, Say Bolt Viscometer Apparatus Tools, Say Bolt Viscometer Apparatus Tool Kits, Say Bolt Viscometer Apparatus Manufacturers, Say Bolt Viscometer Apparatus Suppliers from India, China, Kenya.



Laboratory instruments manufacturers India