

**Product Name :**  
Line-Code Decoder

**Product Code :**  
LIM-CAT-L0017-00006



**Description :**

Line-Code Decoder

**Technical Specification :**

LINE-CODE DECODER  
(PCM WAVEFORMS & DUOBINARY SIGNALING)

Technical Data:

Inputs:

B.CLK TTL level bit clock, synchronized to the input data;  $f_{max} > 100\text{kHz}$

Encoded signal inputs see ENCODER module section of this manual for definitions

Outputs:

DATA decoded TTL level data

STROBE TTL level signal

Line-Code Decoder, Line-Code Decoder Bulk Suppliers, Line-Code Decoder Tools, Line-Code Decoder Line-Code Decoders, Line-Code Decoder Manufacturers, Line-Code Decoder Suppliers from India, China, Kenya.

---

# Laboratory instruments manufacturers India

