

Product Name :
Noise Generator

Product Code :
LIM-CAT-L0017-00007



Description :

Noise Generator

Technical Specification :

Noise Generator

A broadband noise source, with a 12-step output amplitude attenuator.

Technical Data:

Bandwidth 1Hz to Level: 1Vrms at "+22dB" position

Attenuator steps: 12 steps, 0dB to +22dB (2dB per step)

Noise Generator, Noise Generator Bulk Suppliers, Noise Generator Tools, Noise Generator Noise Generators, Noise Generator Manufacturers, Noise Generator Suppliers from India, China, Kenya.