

Product Name :
Fiber Optic Receiver

Product Code :
LIM-CAT-L0017-00025



Description :

Fiber Optic Receiver

Technical Specification :

Fiber Optic Receiver

Technical Data:

Fiber Optic Device: high speed, low noise PIN photo diode, 800nm peak spectral input

Fiber Optic Connector System: single way

Fiber Optic Cable: 1mm, polymer, single core fiber optic cable, sheathed in polyethylene

Output: TTL level digital signal, and, standard level analog signal, switch selectable

Output Frequency Range: DC to >1MHz

Fiber Optic Receiver, Fiber Optic Receiver Bulk Suppliers, Fiber Optic Receiver Tools, Fiber Optic Receiver
Fiber Optic Receivers, Fiber Optic Receiver Manufacturers, Fiber Optic Receiver Suppliers from India, China,
Kenya.

