


Product Name :
Laplace Module

Product Code :
LIM-CAT-L0017-00029



Description :

Laplace Module

Technical Specification :

Laplace Module

Technical Data:

INTEGRATORS: 3-analog-level input integrators with switch adjustable time constants

INTEGRATION RATE: a DIP switch, common to all three integrators, is used to adjust the integration rates to 34V/ms, 24V/ms and 1V/us.

Laplace Module, Laplace Module Bulk Suppliers, Laplace Module Tools, Laplace Module Laplace Modules, Laplace Module Manufacturers, Laplace Module Suppliers from India, China, Kenya.