

Product Name :
Fermented Broth Storage Tank 1500l

Product Code :
LIM-CAT-L0021-00011



Description :

Fermented Broth Storage Tank 1500l

Technical Specification :

Fermented Broth Storage Tank 1500l

Fermented broth storage tank meant to receive the fermented broth after separation of the yeast using a centrifuge separator prior to the distillation process.

It meet the following general specifications:

Material: Stainless Steel

Height: ±2500mm.

Strict leakage and tank pressure tests.

Strictly gasketed sealed tank - no oxygen permeation during storage.

Tank preferably fed from the top and solution to be distilled pumped from the base.

A pumping system with all its components for solution transfer to the distillation pot together with all its spares and accessories.

Easy to operate CIP (Clean-In-Place) system.

Fermented Broth Storage Tank, Fermented Broth Storage Tank Bulk Suppliers, Fermented Broth Storage Tank Tools, Fermented Broth Storage Tank Fermented Broth Storage Tanks, Fermented Broth Storage Tank Manufacturers, Fermented Broth Storage Tank Suppliers from India, China, Kenya.



Laboratory instruments manufacturers India