

**Product Name :**  
Soil Sieve

**Product Code :**  
LIM-CAT-L0042-00010



**Description :**

Soil Sieve

**Technical Specification :**

Soil sieve  
Material: Stainless Steel  
Set of 4  
With mesh sizes 20mm, 2 mm, 0.2 mm, 0.02 mm (gauge set).

Soil Sieve, Soil Sieve Bulk Suppliers, Soil Sieve Tools, Soil Sieve Soil Sieves, Soil Sieve Manufacturers, Soil Sieve Suppliers from India, China, Kenya.