

**Product Name :**  
Light Microscopes

**Product Code :**  
LIM-CAT-L0042-00015



**Description :**

Light Microscopes

**Technical Specification :**

Light Microscopes

Binocular, lens of 4x, 10x, 40x and 100x.

Illumination source : LED (knob with left or right control).

Binocular Microscope With wide field WF10X eyepiece, 4X, 10X, 40X & 100X color coded objectives.

Perfect for various medical, educational, industrial, laboratory and research fields.

With 360 degree Siedentopf type head. Mechanical X-Y axis stage, 120 x 120 mm.

Coaxial coarse and fine adjustment knobs on both sides with 0.002 mm graduations.

LED illumination, comes with dust cover and styrofoam case.

Light Microscopes, Light Microscopes Bulk Suppliers, Light Microscopes Tools, Light Microscopes Light Microscopess, Light Microscopes Manufacturers, Light Microscopes Suppliers from India, China, Kenya.

