

**Product Name :**  
Dissecting Microscopes

**Product Code :**  
LIM-CAT-L0042-00016



**Description :**

Dissecting Microscopes

**Technical Specification :**

Dissecting Microscopes  
20x / 60x dual magnification  
Flexible incident (top) and transmitted (bottom) lights  
Sturdy arm stand with forward widefield optical glass binocular head and all metal framework  
1.3MP USB 2.0 digital camera for photos and videos on Windows XP/Vista/7/8/10 and Mac OSX (with 3rd-party software support)  
User-friendly software (Windows only) offers processing and advanced features including Stitching, EDF, Video Recording and Measurement functions  
45 degree inclined binocular head with locked-in eyepieces.  
Sharp stereo erect images over a wide field view.  
High quality original true color.  
Ocular-tube with adjustable diopter.  
Accurate alignment ensuring a fatigue-free long hour viewing.  
Rugged metal frame and stand.  
Full optical glass elements  
Precise ground glass lens  
All metal mechanical parts  
Adjustable interpupillary distance"  
Built-in both incident and transmitted illumination Systems.

---

LED light.

Approval electric apparatus.

White/black base plate and rubber eye-guards included.

Dissecting Microscopes, Dissecting Microscopes Bulk Suppliers, Dissecting Microscopes Tools, Dissecting Microscopes Dissecting Microscopess, Dissecting Microscopes Manufacturers, Dissecting Microscopes Suppliers from India, China, Kenya.



## Laboratory instruments manufacturers India