

Product Name :
Thermometer

Product Code :
LIM-CAT-L0042-00018



Description :

Thermometer

Technical Specification :

Thermometer
Digital Veterinary Thermometer
For measurement of body temperature of animals
Measuring range temperature: 32...42.9°C
Material: Plastic
Batteries included: Yes
Batteries: 1 x 1.5V LR41
Dimensions: (L) 188 x (W) 26 x (H) 11mm
Weight: 23g.

Thermometer, Thermometer Bulk Suppliers, Thermometer Tools, Thermometer Thermometers, Thermometer Manufacturers, Thermometer Suppliers from India, China, Kenya.

