

Product Name :
Weighing Scales

Product Code :
LIM-CAT-L0042-00019



Description :

Weighing Scales

Technical Specification :

Weighing Scales
Digital 0.01gm-2200gms
18 measuring units including g, kg, oz, and lb
LCD display with backlight and second line
Removable stainless steel weighing pan
USB and RS-232 ports
AC adapter
Auto Shutt Off.
Auto Outside Calibration.
Auto Zero Tracking
Auto Back Light.
Low Battery Indication.
Large Square Weighing Pan.
Net weight Stable indication
Large LCD Display
Dimensions: (Width x Depth x Height) 209x321x98mm.

Weighing Scales, Weighing Scales Bulk Suppliers, Weighing Scales Tools, Weighing Scales Weighing Scales,
Weighing Scales Manufacturers, Weighing Scales Suppliers from India, China, Kenya.



Laboratory instruments manufacturers India