

Product Name :
Eggs Incubator

Product Code :
LIM-CAT-L0042-00020



Description :

Eggs Incubator

Technical Specification :

Eggs incubator

Automatic Egg Turning - Every 2 hours, built-in fan to circulate the air for even temperature and humidity, easy use and clean.

Digital Temperature Control - This egg incubator allows for temperature settings and easy operation which can show important information about temperature, humidity, hatching day and egg turning time.

Cabinet model for egg incubation:

Equipped with 3 layers and movable spacers for different egg sizes.

Each layer can be operated with a separate motor, includes viewing window, interior light and a digital thermometer,

230 volts operation, tough and durable, each tray requires an actuator (not included)

Max. egg output: chicken: 162, goose: 45, duck: 120, turkey: 96 pheasant / guinea fowl, partridge: 252.

Eggs Incubator, Eggs Incubator Bulk Suppliers, Eggs Incubator Tools, Eggs Incubator Eggs Incubators, Eggs Incubator Manufacturers, Eggs Incubator Suppliers from India, China, Kenya.

Laboratory instruments manufacturers India

