

**Product Name :**  
Identification Kit ,e.g Tattooing Kit

**Product Code :**  
LIM-CAT-L0042-00042



**Description :**

Identification Kit ,e.g Tattooing Kit

**Technical Specification :**

Identification kit ,e.g tattooing kit  
Manual plier type model, complete with set of digits A-Z and ink bottles.  
Use on pigs, sheep, goats, cattle and livestock  
Holds up to 4 letters or numbers at a time  
Permanently tattoos the inside of an animal's ear  
High-quality, sandblasted aluminium pliers: lightweight.

Identification Kit ,e.g Tattooing Kit, Identification Kit ,e.g Tattooing Kit Bulk Suppliers, Identification Kit ,e.g Tattooing Kit Tools, Identification Kit ,e.g Tattooing Kit Identification Kit ,e.g Tattooing Kits, Identification Kit ,e.g Tattooing Kit Manufacturers, Identification Kit ,e.g Tattooing Kit Suppliers from India, China, Kenya.

