

**Product Name :**  
Dehorning Rod

**Product Code :**  
LIM-CAT-L0042-00046



**Description :**

Dehorning Rod

**Technical Specification :**

Dehorning Rod

Hot iron type, electric 24V/250W.

Avoid injured by fighting between cows, especially in the breast area is hurt by the follow cow.

Easier to manage , the cow bed and tent area is lesser. Especially the scatter breeding or groups of cattle breeding,

To deal with the Angle of cows is conducted in the early age, often in a cow was born in a week, at the latest not more than two months. So in cows early age, the cow is easy to control, less bleeding, less pain, less susceptible to bacterial infections.

Dehorning Rod, Dehorning Rod Bulk Suppliers, Dehorning Rod Tools, Dehorning Rod Dehorning Rods, Dehorning Rod Manufacturers, Dehorning Rod Suppliers from India, China, Kenya.

