

**Product Name :**  
Multi Temperature Meter

**Product Code :**  
LIM-CAT-L0009-00003



**Description :**

Multi Temperature Meter

**Technical Specification :**

Multi temperature meter

Digital thermometer:

This high quality thermometer is suitable for use with any K-type thermocouple. Each thermometer comes with a 9V battery and spare air probe. It also has following features and specifications:

°C and °F switchable

Measurement range: -50°C – 1300°C (-50°F – 1999°F)

Operating temperature and humidity: 0°C to 50°C (32°F to 122°F), 0-80% RH

Storage temperature and humidity: -10°C to 60°C (14°F to 140°F), 0-70% RH

Data Hold function

T1 & T2 dual points display

T1 – T2 differential measurement

Multi Temperature Meter, Multi Temperature Meter Bulk Suppliers, Multi Temperature Meter Tools, Multi Temperature Meter Multi Temperature Meters, Multi Temperature Meter Manufacturers, Multi Temperature Meter Suppliers from India, China, Kenya.

---

# Laboratory instruments manufacturers India

