

Product Name :
Sieve Shakers Motorized

Product Code :
LIM-CAT-L0044-00002



Description :

Sieve Shakers Motorized

Technical Specification :

Sieve Shakers, Motorized

A digital timer adjustable from 0-99 minutes incorporated as an integral part of the equipment.

The Sieve Shaker is capable of carrying upto 8 sieves.

Driven by a 1/2 HP geared motor Operating on 220 V, single phase AC.

Measuring: 20 µm -25mm

Suitable for Sieving: Dry

Vibrating Frequency: 1400 r/min

Vibrating Amplitude: 1-4mm

Motor: Speed Regulator

Clamping Devices: Supporting clamps & Rods, standard", "comfort", for dry sieving

Housing: Stainless Steel / Powder Coated

Sieve Size Diameter: 200 mm

Quantity of Sieve Capacity: 8 Sieves (200mm x 50mm h)

Laboratory instruments manufacturers India

