

**Product Name :**  
Sieves

**Product Code :**  
LIM-CAT-L0044-00003



**Description :**

Sieves

**Technical Specification :**

Sieves

Sieves are fabricated out of spun brass, without joint frame, double folded bottom, having beading at top, tightly fitting with brass wire mesh.

Brass Sieves: 20 mm, 12.5 mm, 10.00 mm, 6.3mm, 4.75mm, 4.00mm, .80mm, 2.36mm, 2.00mm, 1.70mm, 1.18mm, 1.0mm, 850,600; 425,300, 212, 150, 125, 63microns, pan & cover.