

Product Name :
Automatic Computerized Consolidation Test Apparatus

Product Code :
LIM-CAT-L0044-00015



Description :

Automatic Computerized Consolidation Test Apparatus

Technical Specification :

Automatic Computerized Consolidation Test Apparatus
Vertical force: 20 kN, Ram travel: 25 mm
Testing speed: 0.00001 mm/min,
Testing speed: 50.00000 mm/min,
Horizontal clearance: 175 mm,
Vertical Clearance: 185 mm, (265 mm with extension columns),
Dimension: 300 x 390 x 600 mm
Weights: 40 kg,
Power: 220-110 V, 50-60 Hz, 1 ph.
Additional Accessories: CRS cell, Extension rods, Pressure transducers,
One water pressure line
Cracked Software

