


Product Name :
Temperature Recorder

Product Code :
LIM-CAT-L0009-00005



Description :

Temperature Recorder

Technical Specification :

Temperature recorder

Specifications:

Ambient operating temperature: 0°C to +49°C (+32°F to +120°F). °C or °F selectable.

Chart rotation speed, selectable: 1,5, 3, 6, 12 and 24 hours, 7 and 31 days.

Chart rotation mode, selectable: continuous or single turn.

Chart speed accuracy: ± 1 %.

Chart diameter: 152 mm (6") (60 assorted charts supplied).

Temperature alarm range: -30°C to +50°C (-20°F to +120°F).

Relative humidity alarm range: 0 % to 100 % RH.

Pressure alarm range: 0 bar to 35 bar (0 psi to 500 psi).

Selectable alarm delay (No delay, 10 min, 30 min, 1 h, 2 h) to prevent nuisance alarms.

Remote alarm connection: normally open contacts 48V AC/DC, 0,1 A resistive load.

220 V AC - 240 V AC - 50 Hz / 60 Hz adapter included.

Power consumption: 3,5 W / hour.

Dimensions: 185 mm x 70 mm x 235 mm (7,25" x 2,75" x 9,25"). Weight: 1,13 kg.

Temperature Recorder, Temperature Recorder Bulk Suppliers, Temperature Recorder Tools, Temperature Recorder Temperature Recorders, Temperature Recorder Manufacturers, Temperature Recorder Suppliers from India, China, Kenya.



Laboratory instruments manufacturers India