

**Product Name :**  
Digital PH Meter

**Product Code :**  
LIM-CAT-L0044-00040



**Description :**

Digital PH Meter

**Technical Specification :**

Digital pH meter

Digital PH Meter is on top of people's priority due to their fantastic quality and high durability. They are fabricated with the help of perfect raw materials and international standard technology.

Features consist of: High stability and accuracy Auto facility Slope control 3 1/2 Digit LED display

Display: 3 1/2 Digit LED,

Range pH: -2.00 to 20.00

mV: 0 to  $\pm 1999$ ,

Temperature °C: -5 to 110°C

Resolution PH: 0.01,

Accuracy: 1: 0.01 pH

