

Product Name :
Gas/Liquid Relative Permeameter

Product Code :
LIM-CAT-L0046-00014



Description :

Gas/Liquid Relative Permeameter

Technical Specification :

Pressure Gauge – 0.00 to 1.00 ATM
Golden Core Holder to hold different sizes of core ½” to 1 ½” through rubber ring.
Digital Temperature indicator attached with core holder
A Flow meter (Rota meter) of capacity 0 to 250 cc/sec with regulator holding arrangement.
Regulator to operate different pressure through a gas pressure line attached with pressure gauge.