

Product Name :
Electronic Balance

Product Code :
LIM-CAT-L0047-00006



Description :

Electronic Balance

Technical Specification :

High-Precision digital laboratory gram scale has a measurement 0.01g accuracy. Digital lab scale can meet the requirements of laboratory, kitchen, food, Scientific, and jewelry measurement.

Scientific Scale: It supports 10 weighing units of g/kg/ct/oz/ozt/lb/dwt%/tol/mm, and supports counting mode(pcs), It can be used in factories, laboratories, schools, enterprises, etc.

Food Scale: Large size LCD display, larger numbers, easy to read, more convenient during cooking in the kitchen.

Digital Scale: Our digital laboratory scale has multiple functions such as tare, counting, and memory.