

Product Name :
Lab Electric Stirrers For Mixing Of Drilling Fluids

Product Code :
LIM-CAT-L0047-00011



Description :

Lab Electric Stirrers For Mixing Of Drilling Fluids

Technical Specification :

The range is preferably suited to continuous duty stir and mix applications.
These stirrers are featured with strong adjustable stand, lesser noise and smooth start & stop.

Features:

Heavy duty AC/ DC Universal or Induction Motors.
Coupling is of bush type or self-centring chuck type.
Electronic or Autovariac speed control controllers are additionally given on solicitation.
Coupling is of shrub sort or self-focusing chuck sort.