

Product Name :
Thermostatic Rectangular Water Bath

Product Code :
LIM-CAT-L0047-00013



Description :

Thermostatic Rectangular Water Bath

Technical Specification :

It is a double walled bath with inner chamber made of S.S. & outer wall made of M.S. Sheet duly powder coated/painted.

The gap between outer & inner walls is fitted with special grade glass wool insulation.

The concentric rings can be removed to accommodate containers of various sizes.

Temp. from ambient to $90^{\circ}\text{C} \pm 2^{\circ}\text{C}$ is to be controlled by thermostat.

Specification:

Voltage: 220 Volt AC

Automation Grade: Automatic

Material: Mild Steel

Gain Control: Thermostatically Controlled

Capacity: 8 Litre

Bit Rate: 50 Hz

Power: 220 volt AC

Type: Manual

Accuracy: 1 degree



Laboratory instruments manufacturers India