

Product Name :
Laboratory Sink with Fume Hood

Product Code :
LIM-CAT-L0047-00018



Description :

Laboratory Sink with Fume Hood

Technical Specification :

Fume hoods are made to shield workers from dangerous fumes produced by various processes. These hoods are frequently mounted over equipment where they can filter out hazardous chemicals before they enter the worker's breathing zone.

Fume hoods can be divided into two categories: fixed and portable. Fixed fume hoods are stationary appliances that are indefinitely mounted above machinery. Mobile units called portable fume hoods can be moved around the workspace.

Design Structure: Aerodynamic, Floor mounted

Inside Fume Hood Interior : Polypropylene (PP) Chemical & Heat Resistant, Fire Retardant, Smooth Finish, Easily Laminate 5mm, SS.

Height of worktop from ground level: 850 mm

Airflow Type : Low Constant Volume (for A.C. environment)

Construction (Exterior) :

Anti-corrosive Powder coated

Pre-treated with 8 tank chemical processes and powder coated with highly chemical resistant epoxy Colors having dry film thickness of 70 to 80 microns.

Construction (Interior) : Chemical (especially strong acid like H₂SO₄ and HNO₃) resistant & Heat Resistant, Fire Retardant, Smooth Finish, Easily cleanable Panels



Laboratory instruments manufacturers India