

Product Name :
Wet Test Meters Various Flow Meters Hydrometers

Product Code :
LIM-CAT-L0048-00003



Description :

Wet Test Meters Various Flow Meters Hydrometers

Technical Specification :

Wet test meters, various flow meters hydrometers

Wet Test Meter is a highly accurate gas meters used when accurate measurement of small volumes and low flow rates are necessary.

The Wet Test Meter are highly precise positive displacement style meters. The meters will maintain original calibration for long periods because of the extensive use of corrosion resistant materials.

Applications:

Laboratories and Research Institutions

Small orifice calibration

Determining fuel air ratios for gas appliances

Practical laboratory tool in refinery pilot plant development and research

Measurement of manufactured, natural and specialty (technical) gases

Laboratory experiments

Chemical Reaction gas volumes

Cooking experiment gas usage

Calibration of meters used for atmospheric emission sample testing

Gas Appliance Meter

Calibration Reference.



Laboratory instruments manufacturers India