



Email :
sales@laboratoryinstrumentsmanufacturers.com

Product Name :
Filter Paper For Extraction 240mm

Product Code :
LIM-CAT-L0049-00005



Description :

Filter Paper For Extraction 240mm

Technical Specification :

Filter Paper For Extraction 240mm
Large nominal particle retention for fast analytical filtration.
High wet-strength
High chemical resistance to strong acid and alkali suspensions

Applications:

Use as an animal foodstuff filter for determining fiber content.
Determination of gelatin content in milk and cream.



Laboratory instruments manufacturers India