

Product Name :
Specific Gravity Basket

Product Code :
LIM-CAT-L0049-00015



Description :

Specific Gravity Basket

Technical Specification :

The basket is used in conjunction with a suitable electronic balance, a specially designed robust frame to support the electronic balance and watertight container.

The Density basket is of 20cm Dia x 20cm high fitted with galvanized wire mesh of 6.3mm with handle for suspending it from the balance.

Made of brass with stainless steel wire mesh 6.3mm size ruggedly constructed, 20cm Dia x 20cm high. Complete with handle.