

Product Name :
Water Bath Constant 18*24

Product Code :
LIM-CAT-L0049-00016



Description :

Water Bath Constant 18*24

Technical Specification :

Water Bath Constant 18*24

Rectangular double walled in construction.

Inner chamber made of stainless steel & outer body made of mild steel duly powder coated.

Temperature range 5°C above ambient to 95°C.

Fitted with immersion type heater, digital temperature controller cum indicator accuracy $\pm 0.5^\circ\text{C}$.

Operating on 230 V AC.

Fitted with circular pump with stirrer to re-circulate water at constant temperature.