

**Product Name :**  
Water Bath Constant 18\*24

**Product Code :**  
LIM-CAT-L0049-00016



**Description :**

Water Bath Constant 18\*24

**Technical Specification :**

Water Bath Constant 18\*24

Rectangular double walled in construction.

Inner chamber made of stainless steel & outer body made of mild steel duly powder coated.

Temperature range 5°C above ambient to 95°C.

Fitted with immersion type heater, digital temperature controller cum indicator accuracy  $\pm 0.5^\circ\text{C}$ .

Operating on 230 V AC.

Fitted with circular pump with stirrer to re-circulate water at constant temperature.