

Product Name :
Sectioned Hermetic Compressor

Product Code :
LIM-CAT-L0012-0001



Description :

Sectioned Hermetic Compressor

Technical Specification :

Sectioned Hermetic Compressor

Single-cylinder compressor according to CECOMAF, refrigeration capacity at 5/55°C: 592W, power consumption at 5/55°C: 321W

Evaporating temperature range: -30...10°C

Displacement: 7,95cm³, oil quantity: 450cm³

Refrigerant quantity: 900g, free gas volume in compressor: 1350cm³"

Dimensions and weight

LxWxH: 300x300x240mm

Weight: 12kg

Sectioned Hermetic Compressor, Sectioned Hermetic Compressor Bulk Suppliers, Sectioned Hermetic Compressor Tools, Sectioned Hermetic Compressor Sectioned Hermetic Compressors, Sectioned Hermetic Compressor Manufacturers, Sectioned Hermetic Compressor Suppliers from India, China, Kenya.

Laboratory instruments manufacturers India

