

Product Name :
Electronic Weigh Balance 5kg

Product Code :
LIM-CAT-L0049-00026



Description :

Electronic Weigh Balance 5kg

Technical Specification :

Electronic Weigh Balance 5kg
Large font LCD display with backlight, clear reading.
Stainless Steel Platter
Gram, LB, OZ, CT. Unit
Simple 5 keys panel for easy operation
Adjustable feet
Overload Protection
Below Balance Weighing Facility
Checking Weighing, Piece counting function
External power supply adapts and battery inside.

