

**Product Name :**  
Horizontal Wind Turbine Apparatus

**Product Code :**  
LIM-CAT-L0055-00002



**Description :**

Horizontal Wind Turbine Apparatus

**Technical Specification :**

Blades material: nylon fiber  
Body material: aluminum alloy  
Generator: three phase AC permanent magnet generator/maglev generator  
Control System: electromagnetic brake  
Lubrication mode: self-lubricating  
Phase: Three phase.  
Power= 1kw.  
Stress way of blade= resistance  
Rotating Shaft = Horizontal