

**Product Name :**  
Supply Bench For Training Panel

**Product Code :**  
LIM-CAT-L0058-00001



**Description :**

Supply Bench For Training Panel

**Technical Specification :**

Supply bench for training panels of HL and ST series and others  
Mobile support made of welded square tube, powder coated  
Lockable cabinet  
Drip tray made of glass-fiber reinforced plastic  
Inputs: cold and hot water, gas, mains power and high voltage supplies  
Outputs: 2x 220V, 2x 400V, 1x gas outlet with gas connector.  
1x hot water, 1 x cold water, heating feed and return

**Outputs:**

2x 220V, 50Hz, 1 phase  
1x 400V, 50Hz, 3 phase, 16A  
1x 400V, 50Hz, 3 phase, 32A  
Gas connection with gas connector  
4 hose connections with quick action couplings:  
Heating feed, heating return cold water, hot water

---

# Laboratory instruments manufacturers India

