

Product Name :
Supply Bench For Training Panel

Product Code :
LIM-CAT-L0058-00001



Description :

Supply Bench For Training Panel

Technical Specification :

Supply bench for training panels of HL and ST series and others
Mobile support made of welded square tube, powder coated
Lockable cabinet
Drip tray made of glass-fiber reinforced plastic
Inputs: cold and hot water, gas, mains power and high voltage supplies
Outputs: 2x 220V, 2x 400V, 1x gas outlet with gas connector.
1x hot water, 1 x cold water, heating feed and return

Outputs:

2x 220V, 50Hz, 1 phase
1x 400V, 50Hz, 3 phase, 16A
1x 400V, 50Hz, 3 phase, 32A
Gas connection with gas connector
4 hose connections with quick action couplings:
Heating feed, heating return cold water, hot water

Laboratory instruments manufacturers India

