


**Product Name :**  
Creep Testing Machine

**Product Code :**  
LIM-CAT-L0059-00002



**Description :**

Creep Testing Machine

**Technical Specification :**

Creep Testing machine Table unit for the investigation of the creep specimens made of lead and PP Specimen cross Experiments above or below room temperature possible with conditioning box.

Table unit for the investigation of the creep behaviour of different materials

Flat specimens

Position measurement

Experiments above or below room temperature possible with conditioning box

Force measurement (0.5-50 KN)

Force accuracy (Class 0.5 or ( ? 0.5%)

Center alignment accuracy (? 8%)

Lower pull rod moving speed, range and driven (1.5 mm/min, 200 mm and servo motor 750 w respectively)

Weight (1400 Kg including dead weight).

