


Product Name :
Common Rail Diesel Injection System

Product Code :
LIM-CAT-L0052-00007



Description :

Common Rail Diesel Injection System

Technical Specification :

This is an actual working model of fuel supply system of a CRDI diesel engine. This model helps the student to understand about the parts & working of CRDI system. It is specially designed for demonstration purposes. The model will be supplied with key card & very interesting literature regarding working

Include Training panel with:

Mechanical model of crankshaft and cylinder, Variable rpm, Variable load, Variable temperature, Minimum of four variable injector modes, Regeneration mode, RPM counter, Measurement options on injector, Measurement options at fuel quantity control valve, Measurement options at rail pressure sensor, Measurement options at rail pressure control valve

CD-ROM with Lab software browser and course detail explanation and animation

Storage case for experiment board.

