

Product Name :
Motorized Cut-Away Model Of Demonstration Chassis With
Four Cylinder Engine, RWD

Product Code :
LIM-CAT-L0052-00014



Description :

Motorized Cut-Away Model Of Demonstration Chassis With Four Cylinder Engine, RWD

Technical Specification :

The engine is cut at different locations in order to display its components like piston, cylinder, crank shaft, cam shaft etc.
Cutaway, an original engine with clutch and transmission, fixed on a mobile rack and be driven slowly by means of an electric motor.
This model easily moved in lab on caster wheel and coupled to 220/230V A.C single phase motor through a reduction drive unit.

The model include:
2.5 L, 4 stroke engine, turbo, 4 in-line cylinders
Air Intake System
Fuel system
Exhaust system
Lubrication system
Transmission system
Suspension system
Braking system
Steering mechanism
Electrical system

Differential :Hypoid
Manual with detail explanation.



Laboratory instruments manufacturers India