

**Product Name :**  
Motorized Cut-Away Model Of Demonstration Chassis With  
Four Cylinder Engine, RWD

**Product Code :**  
LIM-CAT-L0052-00014



**Description :**

Motorized Cut-Away Model Of Demonstration Chassis With Four Cylinder Engine, RWD

**Technical Specification :**

The engine is cut at different locations in order to display its components like piston, cylinder, crank shaft, cam shaft etc.  
Cutaway, an original engine with clutch and transmission, fixed on a mobile rack and be driven slowly by means of an electric motor.  
This model easily moved in lab on caster wheel and coupled to 220/230V A.C single phase motor through a reduction drive unit.

The model include:  
2.5 L, 4 stroke engine, turbo, 4 in-line cylinders  
Air Intake System  
Fuel system  
Exhaust system  
Lubrication system  
Transmission system  
Suspension system  
Braking system  
Steering mechanism  
Electrical system

---

Differential :Hypoid  
Manual with detail explanation.



## Laboratory instruments manufacturers India