

**Product Name :**  
Oxygen Cylinder

**Product Code :**  
LIM-CAT-L0054-00004



**Description :**

Oxygen Cylinder

**Technical Specification :**

Type of cylinder MM  
Oxygen capacity: 122 CU - 3455 LITER  
Service Pressure: 2216 PSI - 153 BAR  
Cylinder Length: 35.75.0IN - 90.4CM  
Cylinder: O.D.8 IN. - 203MM  
Cylinder Weight: 38.6 lbs-17.5KG