

Product Name :
Fitting Loss Training Panel

Product Code :
LIM-CAT-L0058-00004



Description :

Fitting Loss Training Panel

Technical Specification :

Training panel for investigating the pressure losses in pipe fittings such as elbows and bends
4 different measuring sections, measured length of each section 2300mm, 10 elbows/bends
Pressure measurement with annular chambers with electronic differential pressure gauge
Flow rate measurement with Rota meter
Hose connections made using quick action coupling
Water feed via pressure reducer
Cold water connection
Differential pressure gauge with bleed
measuring range: 0...2000mbar
supply: 9V, battery-operated
Rota meter: measuring range 150...1600ltr/h
Outlet pressure at pressure reducer: 0.5...2bar
Pipe sections: measured length: 2300mm
Pipe section 1: steel, 1/2", 90° elbow
Pipe section 2: steel, 1/2", 90° bend
Pipe section 3: copper 18x1mm, 90° elbow
Pipe section 4: copper 18x1mm, 90° bend



Laboratory instruments manufacturers India