



Email :
sales@laboratoryinstrumentsmanufacturers.com

Product Name :
Decision Maker

Product Code :
LIM-CAT-L0017-00003



Description :

Decision Maker

Technical Specification :

DECISION MAKER

Technical Data:

Digital waveform inputs two, IN1 and IN2

Digital waveform outputs two, OUT1 and OUT2

Input/Output levels depend upon the waveform format selected:

TTL: +5V, 0V

Unipolar: +2V, 0V

Bipolar: +/-2V

Waveform formats supported:

NRZ-TTL, NRZ-L, NRZ-M, UNI-RZ, BIPOLAR-RZ, RZ-AMI,

Bit Clock input: B.CLK TTL level, nominally 2kHz;

Bit Clock output: B.CLK synchronized to the Output waveform;

Decision point span: > 90% of bit width

DECISION POINT control selection Internal or External, by switch SW2

DECISION POINT control continuous, by front panel knob (INT), or, 0 to 5V DC EXTERNAL input signal (EXT) at Vin

Z-MODULATION level three modes available, with variable level control

Z-MODULATION pulse width 2uS type.

Decision Maker, Decision Maker Bulk Suppliers, Decision Maker Tools, Decision Maker Decision Makers, Decision Maker Manufacturers, Decision Maker Suppliers from India, China, Kenya.



Laboratory instruments manufacturers India