

**Product Name :**  
System Workbench Beech and Steel 1400 x 800 with  
Cabinet

**Product Code :**  
LIM-CAT-L0072-00001



**Description :**

System Workbench Beech and Steel 1400 x 800 with Cabinet

**Technical Specification :**

System Workbench beech/steel 1400 x 800 with Cabinet

Table-top from 40 mm thick multiple-laminated beech veneer, sanded and thinly varnished.

W = 1400 mm, D = 800 mm, H = 750 mm

Steel support, steel cabinet with drawer and compartment with lock and key all metal parts powder coated

Extremely durable and rigid.