

**Product Name :**  
Error Counting Utilities

**Product Code :**  
LIM-CAT-L0017-00004



**Description :**

Error Counting Utilities

**Technical Specification :**

**ERROR COUNTING UTILITIES**

Technical Data:

**EXCLUSIVE-OR GATE**

Inputs A & B TTL level

Output continuous X-OR result or gated with HI time of the input CLK

CLOCK input TTL level,  $f_{max} > 40\text{kHz}$

**MONOSTABLE**

GATE active level DIP switch selectable, active HI or active LO

GATE time:

Normal mode 103, 104, 105, 106

Extended mode normal mode x2, x4 or x8

Expanded mode same as normal or extended modes but divides the pulse count selected by 12

Error Counting Utilities, Error Counting Utilities Bulk Suppliers, Error Counting Utilities Tools, Error Counting Utilities Error Counting Utilities, Error Counting Utilities Manufacturers, Error Counting Utilities Suppliers from India, China, Kenya.



# Laboratory instruments manufacturers India