

**Product Name :**  
Laboratory Microwave Oven

**Product Code :**  
LIM-CAT-L0082-00008



**Description :**

Laboratory Microwave Oven

**Technical Specification :**

Oven is fitted with heating elements on all three sides, ensuring uniform temperature distribution throughout the working chamber

Sealed and enclosed heaters to working chamber for minimum heat loss

Air heaters placed in the air path created at the back for uniform heating and uniform heat distribution Temp. Range.

Special heavy-duty Stainless-steel lock and door hinges with spring and roller mechanism

Stainless steel crimped wire mesh Shelves to withstand heavy load and to offer minimum resistance to air circulation for better heat distribution

Provision to add more shelves, as per requirement of test - Adjustable Shelves to place test samples or Glass wares of varying dimensions

Adjustable air ventilation mounted at top

Operation on 220/230 Volt, 50 Hz, Single phase, AC Supply

With built-in forced air circulation arrangement for uniform heat distribution inside the working chamber

Complete with pilot indicating lamps, 2 meter cord wire and power plug.

---

# Laboratory instruments manufacturers India

