

**Product Name :**  
Electrical Conductivity Meter

**Product Code :**  
LIM-CAT-L0083-00006



**Description :**

Electrical Conductivity Meter

**Technical Specification :**

Digital Conductivity Meter are reliable and accurate test instruments for measurement of Conductivity of aqueous solutions. They measure Conductivity and TDS in five ranges.

The Resolution is  $0.1 \mu\text{S}/\text{Cm}$  and  $0.1 \text{ ppm}$  repetitively, in the lowest range.

The cell constant is flashed on a digital display & can be adjusted from the front panel.

The instruments have the Check facility to calibrate the instrument. These instruments are ideal for monitoring salt contents in water and soluble salts in soils.

They are extremely useful instruments for agriculture and soil analysis labs, swimming pools, fertilizer plants, petroleum refineries, textile plants etc.