

**Product Name :**  
Gas and Bomb Calorimeter

**Product Code :**  
LIM-CAT-L0083-00007



**Description :**

Gas and Bomb Calorimeter

**Technical Specification :**

It is used for determination of calorific value of liquid & Solid fuels.

The Instrument consists S.S. Body made out of corrosion resisting stainless steel Rod cap. 300 ml, Calorimeter vessel, outer Jacket with insulated Body, offset stirrer, S.S. Crucible, pellet press, ignition wire, pressure gauge, with flexible pipes. with electronic firing unit.

Resolution 0.01°C.