

Product Name :
Alternative Energy System Trainer

Product Code :
LIM-CAT-L0084-00006



Description :

Alternative Energy System Trainer

Technical Specification :

NISE/DGET curriculum-based laboratory experiments supported.

Useful for laboratory experimental learning by students in renewable energy basics, energy conservation, charge controller, storage system etc.

Optional single phase (Stand alone or Grid tied) inverter to demonstrate power export using bidirectional multifunction meter.

Closed loop temperature control using peltier module to study temperature effect on solar cell characteristics.