

**Product Name :**  
Fluoride Solution Muffle Furnace

**Product Code :**  
LIM-CAT-L0086-0002



**Description :**

Fluoride Solution Muffle Furnace

**Technical Specification :**

Light weight .run around the chamber & controlled by or digital temperature indicator cum controller, the outer body is made of thick mild steel duly painted. The apparatus is complete with pyrometer. Thermo-couple, silver fuse, two neon lights with wire and plug to work on, (with silicone fibre insulation).

Outer and intermittent walls are made of thick P.C.R.C. sheet, duly powder coated.

Very suitable for various purposes like gravimetric analysis for organic substance and measurement of quantity for coal and lime.

With the Micro Processor Based PID Temp. Controller cum Indicator.

By the highly efficient heater installed in the chamber, it heats the materials uniformly.

With this excellent feature, it surely reduces time of temperature climb.

By using the high density ceramic wool as a insulator, it maintains consistent temperature and minimal heat dissipation.

