

Product Name :
Gas Chromatography

Product Code :
LIM-CAT-L0086-00004



Description :

Gas Chromatography

Technical Specification :

The GC System is now ready to solve your trace analysis needs. This system utilizes the new Barrier Discharge Ionization Detector technology coupled with a GC capillary gas chromatograph to create a GC system that makes it possible to reveal trace components that are difficult to see by other GC detectors.

High Sensitivity

Detection Sensitivity over 100x Higher than TCD, 2x Higher than FID

Novel Universal Detector

Single Detector Approach for Your Complex Analyses

Long-Term Stability

Long-Term Stability with New Discharge Design